

CuroCell UNO®



A new generation of fully automated
alternating pressure mattress replacement systems

CuroCell UNO®

The CuroCell UNO® is used as an aid in the prevention and treatment of pressure ulcers up to and including category 4⁽¹⁾. It stands for innovation and modern design and together with quality engineering it lays the foundation for our new generation of fully automated premium mattress systems.

It is equipped with a large touch screen panel with a user-friendly interface. The fully automated pump offers automatic adjustment, providing optimal pressure relief and comfort for each individual user while aiding simple use.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene⁽²⁾, quality and safety. All covers are manufactured from modern, recyclable materials and are vapour permeable⁽³⁾, which reduce the risk of skin maceration.

Facts

Type: Alternating pressure mattress replacement system

Height: 20 cm

Function: Prevention and treatment up to and including category 4⁽¹⁾

Recommended user weight: 0-250 kg

Warranty: 2 years

Cleaning of cover: Wipe with cleaning agent and/or disinfectants

Machine wash max 95 °C, tumble drying



Features Pump

- Fully automated pump that adapts individually
- Large colour touch display with user-friendly interface
- 3 modes: Alternating, Pulsating and Static
- Pack&Go™ function
- Silent running and vibration free
- Battery backup (optional)
- Soft friction hangers with extra soft tip
- Notification function for loss of air and or electricity supply

Features Mattress

- Cell on Cell construction with safety mattress
- Heel Function
- Static head cells
- CPR-function for rapid deflation
- Smart cable management
- PVC free
- Unique design for easy handling and usage
- Chlorine resistant (10 %) multi stretch cover with welded seams
- Well protected zipper

Sizes

80/85/90/105/120 x 200/210/220 x 20 cm

If you need any other sizes please contact your distributor

Optional

- Durable bag
- Sheet holder

References

(1) National Pressure Ulcer Advisory Panel, European Pressure Ulcer Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers: Quick Reference Guide. Emily Haesler (Ed.). Cambridge Media: Osborne Park, Western Australia; 2014.

(2) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126, SS 876 00 19 "Health care textiles – Bacteria penetration – Wet".

(3) SS-EN ISO 15496:2004, DIN 53.122.



Always read the instruction for use prior to usage of the product.

Care of Sweden and CuroCell UNO are trademarks of Care of Sweden AB.
© Care of Sweden AB, 2015. All rights reserved. Last update: 22/01/15.



Scan QR code to watch presentation video.



Mattress selector

- Guidelines when choosing a product
See how to choose the right mattress at:
www.careofsweden.com/mattress-selector

 care of sweden
Supporting Life

Care of Sweden AB, Box 146, SE-514 23 Tranemo, Sweden
+46 (0)771 106 600 | info@careofsweden.se | www.careofsweden.com